





可拋式金電極

Disposable Gold Electrode




Performance

AUSE100



AUTE100AgCl



Specification

Type	Appearance	Model	Description (Diameter/Area)
Single Electrode		AUSE100	WE : 5mm / 0.196cm ²
Three Electrode		AUTE100	WE : Au (3mm / 0.071cm ²) CE : Au (0.050cm ²) RE : Au (0.010cm ²)
		AUTE100AgCl	WE : Au (3mm / 0.071cm ²) CE : Au (0.050cm ²) RE : AgCl (0.010cm ²)
		AUTE200	WE : Au (2x1mm / 0.020cm ²) CE : Au (0.040cm ²) RE : Au (0.020cm ²)
		AUTE200AgCl	WE : Au (2x1mm / 0.020cm ²) CE : Au (0.040cm ²) RE : AgCl (0.020cm ²)